

**MT1
Mouse Tray**



MAGNATEX

“ORDINARY PEOPLE
DOING EXTRAORDINARY
THINGS”



MAGNATEX, INC.
2520 Ridgemar Court
Louisville, KY 40299
Phone: 502.493.0558
Fax: 502.493.9862
www.magnatexinc.com



MAGNA *TEX* introduces the MT1; Mouse Tray.


The MT1 Mouse Tray can be installed on the Data Acquisition Computer Stand (PWS-02), the Preload Assist Workstation SPA Unloader (PWS-04), and on the Preload Assist Workstation SPA Computer and Monitor Stand (PWS-01A).

The MT1 Mouse Tray easily rotates underneath the keyboard tray for convenient storage.

The MT1 Mouse Tray can be retrofit to any MAGNA *TEX* PWS Keyboard Tray (PWS-01A-004) by drilling a 1/4" hole. Installation is easy with bolt or fastener of your choice.



Contact a MAGNA *TEX* representative today for a complete PARTS, Pre-load assist Units, Conveyors a



For ordering please send an e-mail to gbrown@magnatexinc.com or call us at: 502.493.0558

MAGNATEX is committed to producing and supplying high quality products and services to all our customers.



MagnaTex Inc. meets applicable federal and state program eligibility requirements for several certification programs including 8(a). Certifications available upon request.

The content of this brochure is protected by copyright. You may not copy, distribute, modify, transmit, re-publish for public or commercial purposes without the prior written consent of MagnaTex, Inc. All names, logos, and trademarks that appear in this brochure are the property of MagnaTex, Inc. or are used by MagnaTex Inc. under license. You may not use any of the forgoing for any purpose without the written consent of MagnaTex, Inc.